

TANBead Nucleic Acid Extraction Kits For Maelstrom 9610			
Type of sample	Cat. No.	Description	Sample
Blood	301188	TANBead OptiPure Blood DNA Auto Plate REF: W61EA46; 96 tests/box; Proteinase K included	1. fresh blood* 2. frozen blood*
	301435	TANBead Dried Blood Spot DNA Auto Plate REF: W61EA46-BS ; 96tests/box ; Proteinase K is included	3. buffy coat* 4. parasite in whole blood, serum or plasma*
	301189	TANBead OptiPure Blood DNA Auto Tube REF: W61ES66; 72 tests/box; Proteinase K included	dry blood spot
	301402	TANBead Blood RNA Auto Plate REF: W621A46, 96 tests/box	whole blood
	301403	TANBead Blood RNA Auto Tube REF: W621S66, 72 tests/box	
Animal Tissue	301192	TANBead Tissue Total DNA Auto Plate REF: W6T2A46; 96 tests/box; Proteinase K included	soft or hard animal cells or tissues
	301193	TANBead Tissue Total DNA Auto Tube REF: W6T2S66; 72 tests/box; Proteinase K is included	
	301404	TANBead Tissue Total RNA Auto Plate REF: W6K2A46; 96 tests/box	animal cells or tissues
	301405	TANBead Tissue Total RNA Auto Tube REF: W6K2S66; 72 tests/box	
Plant tissue	301194	TANBead Plant DNA Auto Plate REF: W613A46; 96 tests/box	Plant tissue (leaves)
	301259	TANBead Plant DNA Auto Tube REF: W613S66; 72 tests/box	
	301379	TANBead Plant DNA Auto Plate REF: W613A46-SE ; 96 tests/box ; Proteinase K included	Plant tissue (leaves and seeds)
	301378	TANBead Plant DNA Auto Tube REF: W613S66-SE ; 72 tests/box ; Proteinase K included	
	301406	TANBead Plant RNA Auto Plate REF: W6K3A46; 96 tests/box	Plant tissues
	301407	TANBead Plant RNA Auto Tube REF: W6K3S66; 72 tests/box	
Bacteria	301198	TANBead Gram Bacteria DNA Auto Plate REF: W61GA46 ; 96 tests/box ; Proteinase K and Lysozyme are included	1. Gram-positive and Gram-negative bacteria culture 2. *MTB complex in sputum, BAL, liquid or solid culture
	301199	TANBead Gram Bacteria DNA Auto Tube REF: W61GS66 ; 72 tests/box ; Proteinase K and Lysozyme are included	
	301299	TANBead Gram Bacteria DNA Auto Plate REF: W61GA46-SE ; 96 tests/box ; Proteinase K and Lysozyme are included	Gram-positive and Gram-negative bacteria culture
	301300	TANBead Gram Bacteria DNA Auto Tube REF: W61GS66-SE ; 72 tests/box ; Proteinase K and Lysozyme are included	
Viral	301224	TANBead OptiPure Viral Auto Plate REF: W665A46; 96 tests/box; Proteinase K included	1. serum 2. sputum 3. BAL 4. oropharyngeal swab
	301209	TANBead OptiPure Viral Auto Tube REF: W665S66; 72 tests/box; Proteinase K included	
	301574	TANBead Virapid Virus Auto Plate REF: W685A46, 96 test/box	1. Swab 2. Saliva 3. Urine
	301575	TANBead Virapid Virus Auto Tube REF: W685S46, 96 test/box	
cfDNA	301377	TANBead OptiPure cfDNA Auto Plate REF: W61CA46; 96 tests/box	500-1000 µl blood, plasma, serum
	301386	TANBead OptiPure cfDNA Auto Tube REF: W61CS66; 72 tests/box	
Stool	301392	TANBead Stool Cell DNA Auto Plate W6SCA46; 96 tests/box	stool of human or other animal species
	301391	TANBead Stool Cell DNA Auto Tube W6SCS66; 72 tests/box	